

Product Datasheet

TCF4 antibody (orb413783)

Description	Rabbit polyclonal antibody to TCF4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the C-terminal of Human TCF4 (500-600 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).
Target	TCF4
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-723 was sufficient for detection of TCF4 in 15 µg of mouse liver cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	71 kDa, 45 kDa
Uniprot ID	P15884
NCBI	NP_001077431.1
Entrez	6925
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.