

Product Datasheet

Wnt5a antibody (orb414051)

Description	Rabbit polyclonal antibody against WNT5A conjugated to RPE
Species/Host	Rabbit
Reactivity	Human
Conjugation	RPE
Tested Applications	IHC, WB
Immunogen	Synthetic peptide of Human Wnt5a (200-300 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).
Target	Wnt5a
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-737 was sufficient for detection of Wnt5a in 10 µg of human HeLa cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	42.3 kDa, 35 kDa
Uniprot ID	P41221
NCBI	NP_001243034.1
Entrez	7474
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.