

Product Datasheet

Alpha Synuclein antibody (pSer129) (orb414125)

Description

Rabbit polyclonal antibody to Alpha

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthetic peptide of Human Alpha Synuclein pSer129 (124-134 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).
Target	Alpha Synuclein pSer129
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	0.5mg/mL
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-742 was sufficient for detection of Alpha Synuclein pSer129 in 15 µg of human brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	100, 75, 45, 15 kDa
Uniprot ID	P37840
NCBI	NP_000336.1
Entrez	6622
Dilution Range	WB (1:1000); ICC/IF (1:500)
Expiration Date	12 months from date of receipt.