

Product Datasheet

**Alpha Synuclein antibody (pSer129) [Out of stock]
(orb414126)**

Description

Rabbit polyclonal antibody against

to ATTO-390

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-390

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

Synthetic peptide of Human Alpha Synuclein pSer129 (124-134 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).

Target

Alpha Synuclein pSer129

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

0.5mg/mL

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-742 was sufficient for detection of Alpha Synuclein pSer129 in 15 µg of human brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

100, 75, 45, 15 kDa

Uniprot ID

P37840

NCBI

NP_000336.1

Entrez

6622

Dilution Range

WB (1:1000); ICC/IF (1:500)

Expiration Date

12 months from date of receipt.