

## Product Datasheet

**HSP70 antibody (Acetyl Lys77) [Out of stock]  
(orb414145)**

**Description**

Rabbit polyclonal antibody against

ATTO-488

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

ATTO-488

**Tested**

ELISA, ICC, IHC, WB

**Applications**
**Immunogen**

Synthetic peptide from the N-terminal of Human HSP70 Acetyl Lys77 (1-100 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).

**Target**

HSP70 Acetyl Lys77

**Preservatives**

PBS pH 7.4, 50% glycerol, 0.09% sodium azide \*Storage buffer may change when conjugated

**Concentration**

0.5mg/mL

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

A 1:1000 dilution of SPC-743 was sufficient for detection of HSP70 Acetyl Lys77 in 10 µg of human HeLa cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

**Clonality**

Polyclonal

**MW**

70 kDa

**Uniprot ID**

**P0DMV9, P0DMV8**

**NCBI**

**NP\_005336.3**

**Entrez**

**3303**

**Dilution Range**

WB (1:1000); ELISA (1:1000); ICC/IF (1:60); IHC (1:50)

**Expiration Date**

12 months from date of receipt.