

Product Datasheet

**HSP70 antibody (Acetyl Lys77) [Out of stock]
(orb414147)**

Description

Rabbit polyclonal antibody against

ATTO-594

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

Synthetic peptide from the N-terminal of Human HSP70 Acetyl Lys77 (1-100 aa), conjugated to Keyhole Limpet Haemocyanin (KLH).

Target

HSP70 Acetyl Lys77

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

0.5mg/mL

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-743 was sufficient for detection of HSP70 Acetyl Lys77 in 10 µg of human HeLa cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

70 kDa

Uniprot ID

P0DMV9, P0DMV8

NCBI

NP_005336.3

Entrez

3303

Dilution Range

WB (1:1000); ELISA (1:1000); ICC/IF (1:60); IHC (1:50)

Expiration Date

12 months from date of receipt.