


Product Datasheet

GRIN2D antibody (orb41417)

Description	GRIN2D antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Glutamate receptor ionotropic, NMDA 2D protein (330-480AA)
Target	GRIN2D
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15399
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG-2 whole ce...