



## Product Datasheet

### IFNLR1 antibody (orb41421)


|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | IFNLR1 antibody                                                   |
| <b>Species/Host</b>        | Rabbit                                                            |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                    |
| <b>Immunogen</b>           | Recombinant Human Interferon lambda receptor 1 protein (21-228AA) |
| <b>Target</b>              | IFNLR1                                                            |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                 |
| <b>Form/Appearance</b>     | Liquid                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.     |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG                                                               |
| <b>Clonality</b>           | Polyclonal                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q8IU57</a>                                            |
| <b>Dilution Range</b>      | WB:1:1000-5000, IHC-P:1:20-200                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |



Immunohistochemical staining of human ad...



Immunohistochemical staining of human te...



Western blot analysis of Mouse kidney ti...