

Product Datasheet

NCR2 antibody (orb41433)

Description	NCR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Natural cytotoxicity triggering receptor 2 protein (22-192AA)
Target	NCR2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95944
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.