



Product Datasheet

SLCO2B1 antibody (orb41453)


Description	SLCO2B1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Solute carrier organic anion transporter family member 2B1 protein (130-210AA)
Target	SLCO2B1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94956
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ad...



Immunohistochemical staining of human te...



Western blot analysis of HepG2 whole cel...