

Product Datasheet

UBA1 antibody (Biotin) (orb414612)

Description	UBA1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 1 protein (730-881AA)
Target	UBA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22314
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.