



Product Datasheet

CD3G antibody (orb414657)

Description	CD3G antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human T-cell surface glycoprotein CD3 gamma chain protein (138-182AA)
Target	CD3G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09693
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of HepG2 cell...



Western blot analysis of MCF-7 whole ce...