

Product Datasheet

CYP2A6 antibody (orb414693)

Description	CYP2A6 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cytochrome P450 2A6 protein (86-267AA)
Target	CYP2A6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11509
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.