

Product Datasheet

DUSP1 antibody (orb414709)

Description	DUSP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Dual specificity protein phosphatase 1 protein (125-236AA)
Target	DUSP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28562
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.