



Product Datasheet

ANKS1B antibody (orb41484)


Description	ANKS1B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ankyrin repeat and sterile alpha motif domain-containing protein 1B protein (999-1248AA)
Target	ANKS1B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z6G8
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(lane 1), H...



Immunohistochemical staining of human br...



Immunohistochemical staining of human te...