

Product Datasheet

MCM3AP antibody (Biotin) (orb414848)

Description	MCM3AP antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Germinal-center associated nuclear protein (292-435AA)
Target	MCM3AP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60318
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.