


Product Datasheet

CHRNA antibody (orb41486)

Description	CHRNA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Acetylcholine receptor subunit gamma protein (328-517AA)
Target	CHRNA
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07510
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(lane 1), H1299(lane 2), H1299(lane 3).