



Product Datasheet

TNFRSF9 antibody (orb415025)

Description	TNFRSF9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 9 protein (24-186AA)
Target	TNFRSF9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07011
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of A549 cells...



Western blot analysis of HepG2 whole ce...