

Product Datasheet

TNFRSF9 antibody (FITC) (orb415027)

Description	TNFRSF9 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 9 protein (24-186AA)
Target	TNFRSF9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07011
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.