

Product Datasheet

CACNA1C antibody (FITC) (orb415079)

Description	CACNA1C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Voltage-dependent L-type calcium channel subunit alpha-1C protein (1755-1971AA)
Target	CACNA1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13936
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.