

Product Datasheet

MARK2 antibody (orb415157)

Description	MARK2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein kinase MARK2 protein (373-594AA)
Target	MARK2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7KZI7
Dilution Range	WB:1:2000-10000
Expiration Date	12 months from date of receipt.