


Product Datasheet

PARP2 antibody (orb415237)

Description	PARP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Poly [ADP-ribose] polymerase 2 protein (14-149AA)
Target	PARP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UGN5
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 whole cell...