

Product Datasheet

TUBGCP4 antibody (orb415249)

Description	TUBGCP4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Gamma-tubulin complex component 4 protein (203-309AA)
Target	TUBGCP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UGJ1
Dilution Range	IHC-P:1:20-200, IF:1:200-500
Expiration Date	12 months from date of receipt.