


Product Datasheet

ARID1A antibody (orb415357)

Description	ARID1A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human AT-rich interactive domain-containing protein 1A protein (655-816AA)
Target	ARID1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14497
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...