

Product Datasheet

ATP7A antibody (orb415365)

Description	ATP7A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Copper-transporting ATPase 1 protein (257-488AA)
Target	ATP7A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04656
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.