


Product Datasheet

MOCS2 antibody (orb41539)


Description	MOCS2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Molybdopterin synthase catalytic subunit protein (1-188AA)
Target	MOCS2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O96007
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549(lane 1), H...



Immunohistochemical staining of human ad...



Immunofluorescence analysis of Hela cell...