

Product Datasheet

FUT3 antibody (orb415473)

Description	FUT3 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Galactoside 3(4)-L-fucosyltransferase protein (83-187AA)
Target	FUT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21217
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.