

Product Datasheet

Tubulin beta 3 antibody (HRP) (orb415594)

Description	Tubulin beta 3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Melanocyte-stimulating hormone receptor protein (1-37AA)
Target	MC1R
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01726
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.