


Product Datasheet

MT-ND5 antibody (orb415629)

Description	MT-ND5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human NADH-ubiquinone oxidoreductase chain 5 protein (503-581AA)
Target	MT-ND5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03915
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cell lysate.