


Product Datasheet

NMBR antibody (orb415653)

Description	NMBR antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Neuromedin-B receptor protein (1-41AA)
Target	NMBR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28336
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...