

## Product Datasheet

### SEMA3A antibody (FITC) (orb415847)

|                        |                                                                  |
|------------------------|------------------------------------------------------------------|
| <b>Description</b>     | SEMA3A antibody (FITC)                                           |
| <b>Species/Host</b>    | Rabbit                                                           |
| <b>Reactivity</b>      | Human                                                            |
| <b>Conjugation</b>     | FITC                                                             |
| <b>Immunogen</b>       | Recombinant Human Semaphorin-3A protein (357-445AA)              |
| <b>Target</b>          | SEMA3A                                                           |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b> | Liquid                                                           |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>            | For research use only                                            |
| <b>Isotype</b>         | IgG                                                              |
| <b>Clonality</b>       | Polyclonal                                                       |
| <b>Uniprot ID</b>      | <a href="#">Q14563</a>                                           |
| <b>Dilution Range</b>  | ELISA                                                            |
| <b>Expiration Date</b> | 12 months from date of receipt.                                  |