



Product Datasheet

BLOC1S6 antibody (orb41586)


Description	BLOC1S6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Biogenesis of lysosome-related organelles complex 1 subunit 6 protein (1-170AA)
Target	BLOC1S6
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UL45
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of A431 whole cel...



Immunoprecipitation analysis of HEK293 W...