

Product Datasheet

GDF7 antibody (HRP) (orb415954)

Description	GDF7 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Growth/differentiation factor 7 protein (322-450AA)
Target	GDF7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z4P5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.