



## Product Datasheet

### CHD1L antibody (orb415969)

|                            |                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>         | CHD1L antibody                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                 |
| <b>Reactivity</b>          | Human                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                           |
| <b>Tested Applications</b> | ELISA, IF, WB                                                                          |
| <b>Immunogen</b>           | Recombinant Human Chromodomain-helicase-DNA-binding protein 1-like protein (704-897AA) |
| <b>Target</b>              | CHD1L                                                                                  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                            |
| <b>Form/Appearance</b>     | Liquid                                                                                 |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                          |
| <b>Note</b>                | For research use only                                                                  |
| <b>Isotype</b>             | IgG                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                             |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q86WJ1">Q86WJ1</a>                              |
| <b>Dilution Range</b>      | WB:1:500-5000, IF:1:200-500                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                        |



Immunofluorescence staining of HepG2 cells



Western blot analysis of HepG2 whole cell lysate