

Product Datasheet

CLEC4M antibody (orb416105)

Description	CLEC4M antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human C-type lectin domain family 4 member M protein (238-382AA)
Target	CLEC4M
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H2X3
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.