

Product Datasheet

NR2E3 antibody (Biotin) (orb416128)

Description	NR2E3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Photoreceptor-specific nuclear receptor protein (112-222AA)
Target	NR2E3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5X4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.