

Product Datasheet

FASTKD2 antibody (orb416165)

Description	FASTKD2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human FAST kinase domain-containing protein 2, mitochondrial protein (117-203AA)
Target	FASTKD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYY8
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.