

Product Datasheet

SSX2IP antibody (FITC) (orb416239)

Description	SSX2IP antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Afadin- and alpha-actinin-binding protein (485-614AA)
Target	SSX2IP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2D8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.