



Product Datasheet

GPR32 antibody (orb416297)


Description	GPR32 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Probable G-protein coupled receptor 32 protein (1-44AA)
Target	GPR32
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075388
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-200
Expiration Date	12 months from date of receipt.



IF analysis of A549 cells lysates using ...



Immunohistochemical staining of human ga...



Immunohistochemical staining of human pl...