

Product Datasheet

NUP98 antibody (HRP) (orb416314)

Description	NUP98 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Nuclear pore complex protein Nup98-Nup96 protein (965-1110AA)
Target	NUP98
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P52948
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.