

Product Datasheet

MYO5A antibody (HRP) (orb416398)

Description	MYO5A antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Unconventional myosin-Va protein (1087-1220AA)
Target	MYO5A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y4I1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.