

Product Datasheet

ALG11 antibody (FITC) (orb416439)

Description	ALG11 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human GDP-Man:Man(3)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase protein (416-480AA)
Target	ALG11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2TAA5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.