

Product Datasheet

RAB3GAP2 antibody (orb416449)

Description	RAB3GAP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Rab3 GTPase-activating protein non-catalytic subunit protein (277-387AA)
Target	RAB3GAP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H2M9
Dilution Range	IHC-P:1:200-500
Expiration Date	12 months from date of receipt.