


Product Datasheet

Acetyl-HIST1H4A (K5) antibody (orb416496)

Description	Acetyl-HIST1H4A (K5) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA
Immunogen	Peptide sequence around site of Acetyl-Lys (5) derived from Human Histone H4
Target	HIST1H4A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	ELISA, ChIP
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...