



Product Datasheet

HIST1H3A (phospho-S10) antibody (orb416537)


Description	HIST1H3A (phospho-S10) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	Peptide sequence around site of Phospho-Ser (10) derived from Human Histone H3.1
Target	HIST1H3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunocytochemistry analysis of HeLa cel...



Western blot analysis of HeLa whole cel...



Immunoprecipitation analysis of HeLa who...