


## Product Datasheet

### Acetyl-HIST1H1A (K21) antibody (orb416558)

|                            |                                                                                 |
|----------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>         | Acetyl-HIST1H1A (K21) antibody                                                  |
| <b>Species/Host</b>        | Rabbit                                                                          |
| <b>Reactivity</b>          | Human                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                    |
| <b>Tested Applications</b> | ELISA, IF                                                                       |
| <b>Immunogen</b>           | Peptide sequence around site of Acetyl-Lys (21) derived from Human Histone H1.1 |
| <b>Target</b>              | HIST1H1A                                                                        |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                     |
| <b>Form/Appearance</b>     | Liquid                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                   |
| <b>Note</b>                | For research use only                                                           |
| <b>Isotype</b>             | IgG                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                      |
| <b>Uniprot ID</b>          | <a href="#">Q02539</a>                                                          |
| <b>Dilution Range</b>      | IF:1:50-200                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                 |



Immunofluorescent analysis of HepG2 (sod...