



Product Datasheet

Acetyl-HIST1H1E (K63) antibody (orb416565)


Description	Acetyl-HIST1H1E (K63) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, IF
Immunogen	Peptide sequence around site of Acetyl-Lys (63) derived from Human Histone H1.4
Target	HIST1H1E
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10412
Dilution Range	ICC:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...



Immunofluorescent
analysis of HeLa
cells...



Immunocytochemistry
analysis of HeLa cel...