



Product Datasheet

Acetyl-HIST1H1C (K96) antibody (orb416598)


Description	Acetyl-HIST1H1C (K96) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, ICC, IF, IP, WB
Immunogen	Peptide sequence around site of Acetyl-Lys (96) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	WB:1:500-5000, ICC:1:20-200, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Immunoprecipitation analysis of Hela who...



Immunofluorescent analysis of Hela cells...