



Product Datasheet

Mono-methyl-HIST1H1C (K186) antibody (orb416622)


Description	Mono-methyl-HIST1H1C (K186) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, WB
Immunogen	Peptide sequence around site of Mono-methyl-Lys (186) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Chromatin
Immunoprecipitation
Hela (4*10...



Western blot
analysis of MCF-7
whole ce...



Immunofluorescence
staining of MCF-7
cel...