



Product Datasheet

HIST1H1C (Ab-186) antibody (orb416624)

Description	HIST1H1C (Ab-186) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Peptide sequence around site of Lys (186) derived from Human Histone H1.2
Target	HIST1H1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403
Dilution Range	WB:1:500-5000, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of 293 whole cell lysate



Western blot analysis of HeLa whole cell lysate